

## Adjustable Power Chair Solution

The easy way to secure a wide range of popular power chairs. Automatic hold-down arm is padded to hold securely while cushioning chair. The 28.5" x 38.5" platform (AL560) and 32.5" x 38.5" platform (AL560XL) adjusts for custom fit. Manual crank backup and license plate mount included. Swing-Away option available.

### Features & Benefits

- One Switch Operation
- Adjustable Platform to Provide Custom Fit
- No Chair Modifications Required
- License Plate Mount Included
- Harmar Blue Touch Technology



Specifications	AL560 Automatic Universal Power Chair Lift	AL560XL Extra Large Universal Power Chair Lift
Accommodates	Power Chairs	Larger Power Chairs
Platform	28.5" x 38.5" adjustable platform	32.5" x 38.5" adjustable platform
Lifting Capacity	350 lbs. (160 kg)	350 lbs. (160 kg)
Installed Weight	89 lbs.	101 lbs.
Securement	Hold-down Arm	Hold-down Arm
Hitch	Class II or greater	AL123X Class III hitch adaptor or greater
Warranty	3-year transferable	3-year transferable